

# Application Information

## Ultracentrifugation

### A Rapid Method for Ribosome Preparation: Part 1— Using the High-Capacity Fixed Angle MLA-80 Rotor in an Optima™ MAX Tabletop Ultracentrifuge

Paul Voelker, Donald K. McRorie, Melvin Dorin  
Beckman Coulter, Inc.

#### Introduction

The separation performance of the Optima™ MAX ultracentrifuge is demonstrated by isolating a bacterial ribosome using two high-capacity rotors. This work is described in two separate application information bulletins:

- Part 1 (this document) describes conditions for pelleting an intact 70 S ribosome using differential centrifugation with a fixed angle MLA-80 rotor and an 8-mL tube.
- Part 2 (document number A-1851A) describes conditions for preparing and separating 30 S and 50 S ribosomal subunits using rate zonal density gradient centrifugation with a swinging bucket MLS-50 rotor and a 5-mL tube.

Run conditions for part 1 included a final pelleting spin of 1 hr 50 min at  $444,000 \times g$  (80,000 rpm) and 5°C. Results show that approximately 350 mg of ribosomes can be isolated using this rotor and tube combination. A pelleting time of less than two hours allows for the processing of a ribosome

preparation, from bacterial cells to final isolation, in a single day. The ribosome was characterized using analytical ultracentrifugation (Optima XL-I from Beckman Coulter) and electron microscopy.

#### Method

Samples of the intact 70 S species of ribosome were prepared from *Escherichia coli* (migula strain) following standard procedures.<sup>(1,2,3)</sup> An example of a step-wise approach used to isolate the ribosome is presented below. All procedures were performed at 5°C unless otherwise indicated.

#### Buffers:

- A: 20 mM MgCl<sub>2</sub> 200 mM NH<sub>4</sub>Cl, 6 mM 2-mercaptoethanol, 0.1 mM Na<sub>2</sub>EDTA, 10 mM Tris-HCl (pH 7.4)
- B: 10 mM MgCl<sub>2</sub> 500 mM NH<sub>4</sub>Cl, 6 mM 2-mercaptoethanol, 0.1 mM Na<sub>2</sub>EDTA, 10 mM Tris-HCl (pH 7.4)
- C: 10 mM MgCl<sub>2</sub> 50 mM NH<sub>4</sub>Cl, 6 mM 2-mercaptoethanol, 0.1 mM Na<sub>2</sub>EDTA, 10 mM Tris-HCl (pH 7.4)


1. *E. coli* (host for plasmid vector pBR322) were grown overnight in 2 L of LB media at 37°C with vigorous shaking (300 rpm).
2. The cells were harvested in a Avanti® J-30I high-speed centrifuge from Beckman Coulter with a JLA-8.1000 rotor using two 1 L bottles at  $10,000 \times g$  for 6 minutes, and the pellets pooled.
3. The bacterial pellet was resuspended in an equal volume of buffer A, then ruptured in a cell disruption bomb using nitrogen at 2,500 psi. To maximize the percentage of ruptured cells, the cells were passed through the bomb twice.
4. A solution of 1  $\mu\text{g/mL}$  of DNase in water was added to the cell homogenate.
5. The cell homogenate was divided into two aliquots and spun in two 50 mL bottles at  $105,000 \times g$  for 4 hours in an Avanti® J-30I centrifuge with a JA-30.50 rotor from Beckman Coulter. The supernatant was discarded.
6. The pellet from each bottle was resuspended in buffer A and split between two 8 mL tubes. A 4-mL solution of 30% sucrose in buffer B was added to each 8-mL Quick-Seal® tube, followed by careful addition of approximately 4 mL of the ribosome suspension to the layer. The material was centrifuged in an MLA-80 rotor for 1 hr 50 min at  $444,000 \times g$  (80,000 rpm). (Note: the tube contained a visible contaminant that banded approximately two-thirds of the way from the top of the tube.)
7. The supernatant was decanted and the clear pellet resuspended in approximately 1 mL of buffer C in a 50-mL bottle and clarified by spinning at  $20,000 \times g$  for 30 minutes in an Avanti J-30I centrifuge with a JA-30.50 rotor. The supernatant was dialyzed against buffer C and stored at 5°C.

The physical properties of the separation were determined with an Optima™ XL-I analytical ultracentrifuge from Beckman Coulter using a sedimentation velocity protocol. Runs were made at 30,000 rpm at 25°C in buffer C. Data were measured using interference optics and size distributions were determined using the time-derivative analysis of the sedimentation velocity profiles.<sup>(4)</sup> Further characterization of the ribosome preparation was accomplished using electron microscopy (EM). To preserve the


ribosome prior to EM analysis, a sample of the pellet was treated with a solution of 2.5% glutaraldehyde (EM grade) in 0.1 M sodium phosphate buffer (pH 7.4) and stored at 5°C. The pellet was post-fixed with 2% osmium tetroxide and stained with 0.5% uranyl acetate and Reynolds lead stain. Ultrathin sections were viewed and photographed on a Zeiss model 902 microscope. Electron microscopy was performed by Dr. Megan Yao at Stanford Research Institute (SRI), Menlo Park, CA.

## Results/Discussion

Analysis of the separation of bacterial ribosome showed the final pellet contained a majority of the 70 S ribosomes, with dissociated 30 S and 50 S subunits present only as a small percentage (Figure 1). The ribosomes provided a ratio of 2.04 for  $A_{260}/A_{280}$ . (A ratio  $> 1.8$  for  $A_{260}/A_{280}$  is an accepted criterion for a pure ribosome preparation).<sup>(5)</sup> Electron microscopy confirmed that the isolate was highly enriched in ribosomes (Figure 2).



**Figure 1.** Isolation of bacterial ribosomes following centrifugation (Table 1). The size distributions are presented graphically as peaks and the preparation is shown to contain a majority of the intact 70 S with the two subunits, 30 S and 50 S, present in smaller relative yields—estimated by measuring the areas under the peaks. Sedimentation coefficients ( $s^*$ ) were measured at 20°C and are corrected for viscosity. Error bars are included on the graph.



**Figure 2.** Electron micrograph of 70 S ribosomes following centrifugation ( $\times 140,000$ ). The shape and size of the isolated ribosomes are consistent with expectations. Bar represents  $0.05 \mu\text{m}$ .

Approximately 43 mg of ribosomes were recovered from an MLA-80 rotor using a single 8-mL tube. Considering the MLA-80 has an 8-tube capacity, a maximum of 350 mg can be separated using this rotor and tube combination. Previously, the largest tube volume available for a tabletop ultracentrifuge was 5.1 mL (used with a TLA-100.4). The increased tube capacity of the MLA-80 rotor allows for a greater throughput than was possible previously. Larger tubes allow more material to be pooled and will reduce losses. Run times of less than 2 hours provide the added advantage of being able to isolate a ribosome from a bacterial mixture in a single day. Results are presented in Table 1.

**Table 1. Process Capacity of a Bacterial Ribosome Using an Optima™ MAX Tabletop Ultracentrifuge**

	<i>MLA-80</i> *
Loading	32 mL homogenate ( $8 \times 8 \text{ mL}$ )†
Recovery	~344 mg
Run Time	1 hr 50 min

\*  $444,000 \times g$  at  $5^\circ\text{C}$ .

† In each tube, 4 mL of homogenate was layered over 4 mL of a 30% sucrose solution.

## References

1. Maniatis, T., Fritsch, E. F., Sambrook, J. *Molecular Cloning*. Cold Spring Harbor, NY: Cold Spring Harbor Laboratory, 1982.
2. Noller, H. F., Jr., and Moldave, K., Editors. *Methods in Enzymology*, vol. 164, "Ribosomes," pp. 428-430. San Diego: Academic Press, Inc., 1988.
3. Fraser, D. Bursting Bacteria by Release of Gas Pressure. *Nature* 167, 33-34 (1951)
4. Stafford, W. F. Boundary Analysis in Sedimentation Transport Experiments: A Procedure in Obtaining Sedimentation Coefficient Distributions Using the Time Derivative of the Concentration Profile. *Anal. Biochem.* 203, 295-301 (1992)
5. G. Spedding, Editor. *Ribosomes and Protein Synthesis—A Practical Approach*, pp. 5-6. Oxford: IRL Press at Oxford University Press, 1990.

All trademarks are the property of their respective owners.



Beckman Coulter, Inc. • 4300 N. Harbor Boulevard, Box 3100 • Fullerton, California 92834-3100

Sales: 1-800-742-2345 • Service: 1-800-551-1150 • Internet: [www.beckmancoulter.com](http://www.beckmancoulter.com) • Telex: 678413 • Fax: 1-800-643-4366

---

**Beckman Coulter Worldwide Bioresearch Division Offices:**

**Africa, Middle East, Eastern Europe (Switzerland)** (41) 22 994 07 07 **Australia** (61) 2 9844-6000 **Austria** (43) 1 729 2164 **Canada** (905) 819-1234  
**China** (8610) 6527 9344-9 **France** (33) 1 43 01 70 00 **Germany** (49) 89 35870-0 **Hong Kong** (852) 2814 7431 **Italy** (39) 2-953921  
**Japan** 03-5352-2820 **Mexico** 525-559-16-35 **Netherlands** (31) 297-230630 **Singapore** (65) 339 3633 **South Africa** (27) 11-805-2014/5  
**Spain** (34) 1358 0051 **Sweden** (46) 898-5320 **Switzerland** (41) 22 994 07 46 **Taiwan** (886) 02 2378-3456 **U.K.** (44) 1494 441181  
**U.S.A.** 1-800-742-2345.

