

# Application Information

## Ultracentrifugation

### A Rapid Method for Ribosome Subunit Isolation: Part 2— Using the High-Capacity Swinging Bucket MLS-50 Rotor in an Optima™ MAX Tabletop Ultracentrifuge

*Paul Voelker and Donald K. McRorie  
Beckman Coulter, Inc.*

#### Introduction

The separation performance of the Optima™ MAX ultracentrifuge is demonstrated by isolating a bacterial ribosome using two high-capacity rotors. This work is described in two separate application information bulletins:

- Part 1 (document number A-1850A) describes conditions for pelleting an intact 70 S ribosome using differential centrifugation with a fixed angle MLA-80 rotor and an 8-mL tube.
- Part 2 (this document) describes conditions for preparing and separating 30 S and 50 S ribosomal subunits using rate zonal density gradient centrifugation with a swinging bucket MLS-50 rotor and a 5-mL tube.

Results for part 2 show approximately 3 mg of dissociated ribosomes are separated as 30 S and 50 S subunits with negligible cross-contamination using this rotor/tube combination. A 5-mL tube replaces the 2.2-mL tube as the largest available size and allows for greater throughput than was possible previ-

ously in a tabletop ultracentrifuge. The subunits were analyzed by analytical ultracentrifugation.

#### Method

The 70 S ribosome was dissociated and separated into the 30 S and 50 S subunits following standard procedures.<sup>(1,2,3)</sup> (Note: preparation of the intact 70 S species is described in Application Bulletin A-1850A.) An example of a step-wise approach used to prepare and isolate ribosomal subunits is presented below. All procedures were performed at 5°C unless otherwise indicated.

#### Buffers


- D: 20 mM MgCl<sub>2</sub> 50 mM NH<sub>4</sub>Cl,  
6 mM 2-mercaptoethanol, 10 mM Tris-HCl  
(pH 7.4)
- E: 20 mM MgCl<sub>2</sub> 100 mM NH<sub>4</sub>Cl,  
6 mM 2-mercaptoethanol, 10 mM Tris-HCl  
(pH 7.4)
- F: 1 mM MgCl<sub>2</sub> 500 mM NH<sub>4</sub>Cl,  
6 mM 2-mercaptoethanol, 10 mM Tris-HCl  
(pH 7.4)

1. A bacterial pellet was resuspended in an equal volume of buffer D, then ruptured in a cell disruption bomb using nitrogen at 2500 psi. To maximize the percentage of ruptured cells, the cells were passed through the bomb twice.
2. A solution of 1  $\mu\text{g}/\text{mL}$  of DNase (in water) was added to the cell homogenate.
3. The cell homogenate was divided into two aliquots and spun in two 50 mL bottles at  $105,000 \times g$  for 4 hours in an Avanti<sup>®</sup> J-30I centrifuge from Beckman Coulter with a JA-30.50 rotor. The supernatant was discarded.
4. The pellet was resuspended in the same buffer, clarified for 30 minutes at  $20,000 \times g$ , and precipitated with ammonium sulfate (49 g of solid ammonium sulfate per 100 mL of ribosome suspension with a concentration of approximately 1 mg/mL) and stored overnight at 4°C.
5. The precipitate was collected by spinning for 30 minutes at  $30,000 \times g$ .
6. The resulting ribosomal pellet was resuspended in a small volume of buffer E, dialyzed against the same buffer, and then against buffer F for dissociation into subunits.
7. The suspension was again clarified by spinning for 30 minutes at  $30,000 \times g$ . An additional pelleting step was required to concentrate the solution of subunits into a smaller volume—the solution was centrifuged in an MLS-50 rotor at  $268,000 \times g$  (50,000 rpm) for 5 hr 15 min at 5°C. The pellet was resuspended in 1 mL of buffer F and stored at 5°C.
8. The individual subunits were isolated by layering the ribosomal suspension (in a 250- $\mu\text{L}$  volume) onto a 5-mL tube containing a linear 25-40% sucrose gradient in buffer F and spinning in an MLS-50 rotor at  $268,000 \times g$  (50,000 rpm) for 5 hr at 5°C.
9. Following centrifugation, the tube bottom was pierced and fractions of approximately 20 drops were collected. Fractions were dialyzed against buffer F and stored at 5°C.

The physical properties of 1) dissociation and 2) separation were determined with an Optima<sup>™</sup> XL-I analytical ultracentrifuge from Beckman Coulter using a sedimentation velocity protocol. Runs were made at 30,000 rpm at 20°C in buffer F. Data were measured using interference optics and the size distributions were determined using the time derivative analysis of the sedimentation velocity profiles.<sup>(4)</sup> Data are presented graphically as peaks. The areas under the peaks are used to estimate the relative yield of each subunit.

## Results/Discussion


The dissociation of bacterial ribosomes (70 S) produced the two subunits (30 S and 50 S) in near quantitative yield (Figure 1). A 2:1 ratio of 50 S to 30 S subunits (measured as weight % by the XL-I) with a  $A_{260}/A_{280}$  ratio of 1.95 provided an estimate of sample purity (a ratio for  $A_{260}/A_{280} > 1.8$  is an acceptable criterion for a pure ribosome preparation).<sup>(5)</sup>



**Figure 1.** Dissociation of bacterial ribosomes (70 S) into the 50 S and 30 S subunits. A 2:1 ratio of 50 S to 30 S subunits (measured as weight % by the Optima<sup>™</sup> XL-I) was estimated by measuring the areas under the peaks. Sedimentation coefficients ( $s^*$ ) were measured in aqueous buffer at 20°C.

The two ribosomal subunits were separated in 5 hours using a 25-40% sucrose density gradient in a 5-mL tube with a swinging bucket MLS-50 rotor at 5°C. (A 10-40% gradient produced band separation near the tube bottom rather than the middle under the same conditions.) A swinging bucket rotor is considered the rotor of choice for rate zonal density gradient centrifugation. At speed, the horizontal orientation of the bucket ensures that a sample moves through the gradient without encountering any disruptive wall effects.

Sample capacity of the 5-mL tube was determined by loading various amounts of a ribosomal mixture to the gradient and measuring the resolution between the two subunits. The sample mixture was quantified by measuring the absorbance at 260 nm (one  $A_{260}$  unit equals approximately 0.26 mg of bacterial pellet—using reported values of 69 pmol and 34.5 pmol for 30 S and 50 S subunits, respectively).<sup>(5)</sup> Optimal sample capacity was found to be approximately 3 mg with a loading volume of



**Figure 2.** Isolation of a ~3 mg mixture of ribosomal subunits by rate zonal density gradient centrifugation. The 30 S (left-hand plot) and 50 S (right-hand plot) subunits were separated with negligible cross-contamination. Sedimentation coefficients ( $s^*$ ) were measured in aqueous buffer at 20°C.

250  $\mu$ L. Analysis of the separation showed the ribosomal mixture was isolated as two distinct subunits with negligible cross-contamination (Figure 2). The two subunits were collected in three successive fractions corresponding to a separation occurring near the middle of the tube. Table 1 shows that the 50 S and 30 S subunits were recovered in quantities of 1.85 mg and 0.85 mg, respectively, with an overall yield of 2.7 mg (86%). A ratio of 2:1 at  $A_{260}$  for 50 S to 30 S and a ratio of about 2 for  $A_{260}/A_{280}$  for each of the subunits were consistent with expectations.

**Table 1. Rate Zonal Separation of Bacterial Ribosomal Subunits Using an Optima™ MAX Tabletop Ultracentrifuge<sup>1</sup>**

	Recovery (mg)	$A_{260}$	$A_{280}$
70 S <sup>2</sup> (starting material)	3.13	0.98	0.5
30 S <sup>3</sup>	1.85	0.33	0.15
50 S <sup>3</sup>	0.85	0.72	0.33
(total)	2.7		

1. The two subunits were separated in 5 hours using a 25-40% sucrose density gradient in a 5-mL tube with an MLS-50 rotor.
2.  $A_{260}$  and  $A_{280}$  measurements correspond to a dilution of 1:100.
3.  $A_{260}$  and  $A_{280}$  measurements correspond to a dilution of 1:10.

The separation was repeated with 6 mg of a ribosomal mixture (also in a 250  $\mu$ L volume) loaded on the gradient. Results (not shown) produced a separation with significant cross-contamination of the two subunits. In addition, the two subunits were distributed over a larger number of tube fractions (four instead of three). A separation with 1.5 mg of ribosomal mixture was comparable in quality to a 3 mg sample load but made use of only half the loading capabilities of the tube.

## Conclusions

Before introduction of the MLS-50 rotor, separation of subcellular organelles (and other biological material) using rate zonal density gradient centrifugation in a swinging bucket rotor was generally performed in a floor model ultracentrifuge because of the availability of larger-capacity rotors. With the introduction of the MLS-50 swinging bucket rotor for the Optima MAX tabletop ultracentrifuge, this type of separation is possible using a smaller instrument. The availability of a 5-mL tube in a tabletop instrument provides additional capacity over a 2.2-mL tube and allows for separations on the same scale as those obtained in rotors in a floor model ultracentrifuge.

## References

1. Maniatis, T., Fritsch, E. F., Sambrook, J.  
*Molecular Cloning*. Cold Spring Harbor, NY:  
Cold Spring Harbor Laboratory, 1982.
2. Noller, H. F., Jr., and Moldave, K., Editors.  
*Methods in Enzymology*, vol. 164, "Ribosomes,"  
pp. 428-430. San Diego: Academic Press, Inc.,  
1988.
3. Fraser, D. Bursting Bacteria by Release of Gas  
Pressure. *Nature* 167, 33-34 (1951)
4. Stafford, W. F. Boundary Analysis in Sedimenta-  
tion Transport Experiments: A Procedure in  
Obtaining Sedimentation Coefficient Distributions  
Using the Time Derivative of the Concentration  
Profile. *Anal. Biochem.* 203, 295-301 (1992)
5. G. Spedding, Editor. *Ribosomes and Protein Syn-  
thesis—A Practical Approach*, pp. 5-6. Oxford:  
IRL Press at Oxford University Press, 1990.

---

*All trademarks are the property of their respective owners.*



Beckman Coulter, Inc. • 4300 N. Harbor Boulevard, Box 3100 • Fullerton, California 92834-3100  
Sales: 1-800-742-2345 • Service: 1-800-551-1150 • Internet: [www.beckmancoulter.com](http://www.beckmancoulter.com) • Telex: 678413 • Fax: 1-800-643-4366

---

**Beckman Coulter Worldwide Bioresearch Division Offices:**

**Africa, Middle East, Eastern Europe (Switzerland)** (41) 22 994 07 07 **Australia** (61) 2 9844-6000 **Austria** (43) 1 729 2164 **Canada** (905) 819-1234  
**China** (8610) 6527 9344-9 **France** (33) 1 43 01 70 00 **Germany** (49) 89 35870-0 **Hong Kong** (852) 2814 7431 **Italy** (39) 2-953921  
**Japan** 03-5352-2820 **Mexico** 525-559-16-35 **Netherlands** (31) 297-230630 **Singapore** (65) 339 3633 **South Africa** (27) 11-805-2014/5  
**Spain** (34) 1358 0051 **Sweden** (46) 898-5320 **Switzerland** (41) 22 994 07 46 **Taiwan** (886) 02 2378-3456 **U.K.** (44) 1494 441181  
**U.S.A.** 1-800-742-2345.