Q&A

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REFERENCES

- 1 Sah S, Chen L, Houghton J, et al. Functional DNA quantification guides accurate next-generation sequencing mutation detection in formalin-fixed, paraffin-embedded tumor biopsies. *Genome Med* 2013;5:77.
- 2 Blow N. Tissue preparation: tissue issues. *Nature* 2007;448:959-63.
- 3 Li P, Conley A, Zhang H, et al. Whole-Transcriptome profiling of formalin-fixed, paraffin-embedded renal cell carcinoma by RNA-seq. *BMC Genomics* 2014;15:1087.
- 4 Do H, Dobrovic A. Dramatic reduction of sequence artefacts from DNA isolated from formalin-fixed cancer biopsies by treatment with uracil- DNA glycosylase. *Oncotarget* 2012;3:546-558.
- 5 Kerick M, Isau M, Timmermann B, et al. Targeted high throughput sequencing in clinical cancer settings: formaldehyde fixed-paraffin embedded (FFPE) tumor tissues, input amount and tumor heterogeneity. *BMC Med Genomics* 2011;4:68.
- 6 Hedegaard J, Thorsen K, Lund MK, et al. Next-generation sequencing of RNA and DNA isolated from paired fresh-frozen and formalin-fixed paraffin-embedded samples of human cancer and normal tissue. *PLoS One* 2014;9:e98187.
- 7 Astolfi A, Urbini M, Indio V, et al. Whole exome sequencing (WES) on formalin-fixed, paraffin-embedded (FFPE) tumor tissue in gastrointestinal stromal tumors (GIST). *BMC Genomics* 2015;16:892.
- 8 Graw S, Meier R, Minn K, et al. Robust gene expression and mutation analyses of RNA-sequencing of formalin-fixed diagnostic tumor samples. *Sci Rep* 2015;5:12335.
- 9 Spencer DH, Sehn JK, Abel HJ, et al. Comparison of clinical targeted next-generation sequence data from formalin-fixed and fresh-frozen tissue specimens. *J Mol Diagn* 2013;15:623-633.
- 10 Norton N, Sun Z, Asmann YW, et al. Gene expression, single nucleotide variant and fusion transcript discovery in archival material from breast tumors. *PLoS One* 2013;8:e81925.
- 11 Munchel S, Koestler DC, Bibikova M. Targeted or whole genome sequencing of formalin fixed tissue samples: Potential applications in cancer genomics. *Oncotarget* 2015;6(28):25943-61.
- 12 Blum F: Der Formaldehyde als Hartungsmittel. *Z wiss Mikr* 1893;10:314.
- 13 Dove A. Hard-core sequencing. Science/AAAS Custom Publishing Office. 2016. Available from: URL: <http://www.sciencemag.org/custom-publishing/technology-features/hard-core-sequencing>.
- 14 Feldman MY. Reactions of nucleic acids and nucleoproteins with formaldehyde. *Prog Nucleic Acid Res Mol Biol* 1973;13:1- 49.
- 15 Sikorsky JA, Primerano DA, Fenger TW, et al. DNA damage reduces Taq DNA polymerase fidelity and PCR amplification efficiency. *Biochem Biophys Res Commun* 2007;355(2):431-437.
- 16 Bruskov VI, Malakhova LV, Masalimov ZK, et al. Heat-induced formation of reactive oxygen species and 8-oxoguanine, a biomarker of damage to DNA. *Nucleic Acids Res* 2002;30(6):1354-1363.
- 17 Pfeifer GP, You YH, Besaratinia A. Mutations induced by ultraviolet light. *Mutat Res* 2005;571(1-2):19-31.
- 18 Wagle N, Berger MF, Davis MJ, et al. High-throughput detection of actionable genomic alterations in clinical tumor samples by targeted, massively parallel sequencing. *Cancer Discov* 2012;2:82-93.
- 19 Yost SE, Smith EN, Schwab RB, et al. Identification of high-confidence somatic mutations in whole genome sequence of formalin-fixed breast cancer specimens. *Nucleic Acids Res* 2012;40(14):e107.
- 20 Wong SQ, Li J, Tan AY, et al. Sequence artefacts in a prospective series of formalin-fixed tumours tested for mutations in hotspot regions by massively parallel
- 21 Masuda N, Ohnishi T, Kawamoto S, et al. Analysis of chemical modification of RNA from formalin-fixed samples and optimization of molecular biology applications for such samples. *Nucleic Acids Res* 1999;27:4436-4443.
- 22 Do H, Wong SQ, Li J, Dobrovic A. Reducing sequence artifacts in amplicon-based massively parallel sequencing of formalin-fixed paraffin-embedded DNA by enzymatic depletion of uracil-containing templates. *Clin Chem* 2013;59:1376 - 83.
- 23 Anaka M, Hudson C, Lo PH, et al. Intratumoral genetic heterogeneity in metastatic melanoma is accompanied by variation in malignant behaviors. *BMC Med Genomics* 2013;6:40.
- 24 Do H, Dobrovic A. Limited copy number-high resolution melting (LCN-HRM) enables the detection and identification by sequencing of low level mutations in cancer biopsies. *Molec Cancer* 2009;8:82.
- 25 Schweiger MR, Kerick M, Timmermann B, et al. Genome-wide massively parallel sequencing of formaldehyde fixed-paraffin embedded (FFPE) tumor tissues for copy-number- and mutation-analysis. *PLoS ONE* 2009;4:e5548.
- 26 Sengüven B, Baris E, Oygur T, et al. Comparison of methods for the extraction of DNA from formalin-fixed, paraffin-embedded archival tissues. *Int J Med Sci* 2014;11(5):494-99.
- 27 Bolognesi C, Forcato C, Buson G, et al. Digital Sorting of Pure Cell Populations Enables Unambiguous Genetic Analysis of Heterogeneous Formalin-Fixed Paraffin-Embedded Tumors by Next Generation Sequencing. *Sci Rep* 2016;6:20944.
- 28 Kokkat TJ, Patel MS, McGarvey D, et al. Archived formalin-fixed paraffin-embedded (FFPE) blocks: A valuable underexploited resource for extraction of DNA, RNA, and protein. *Biopreserv and Biobank* 2013;11(2):101-06.
- 29 Corless CL, Spellman PT. Tackling formalin-fixed, paraffin-embedded tumor tissue with next-generation sequencing. *Cancer Discov* 2012;2:23-24.
- 30 Ludyga N, Grünwald B, Azimzadeh O, et al. Nucleic acids from long-term preserved FFPE tissues are suitable for downstream analyses. *Virchows Arch* 2012;460(2):131-40.
- 31 Gilbert MT, Haselkorn T, Bunce M, et al. The isolation of nucleic acids from fixed, paraffin-embedded tissues-which methods are useful when? *PLoS One* 2007;2(6):e537.
- 32 Wang M, Escudero-Ibarz L, Moody S, et al. Formalin-Fixed, Paraffin-Embedded Tissues by Fluidigm Multiplex PCR and Illumina Sequencing. *J Mol Diagn* 2015;17:521e532.
- 33 Ausch C, Buxhofer-Ausch V, Oberkanins C, et al. Sensitive detection of KRAS mutations in archived formalin-fixed paraffin-embedded tissue using mutant-enriched PCR and reverse-hybridization. *J Mol Diagn* 2009;11:508-513.
- 34 Solassol J, Ramos J, Crapez E, et al. KRAS Mutation Detection in Paired Frozen and Formalin-Fixed Paraffin-Embedded (FFPE) Colorectal Cancer Tissues. *Internat J Mol Sci* 2011;12:3191-3204.
- 35 Beltran H, Yelensky R, Frampton GM, et al. Targeted Next-generation Sequencing of Advanced Prostate Cancer Identifies Potential Therapeutic Targets and Disease Heterogeneity. *Eur Urol* 2013;63:920-26.
- 36 Van Allen EM, Wagle N, Stojanov P, et al. Whole-exome sequencing and clinical interpretation of formalin-fixed, paraffin-embedded tumor samples to guide precision cancer medicine. *Nat Med* 2014;20:682-8.
- 37 Wood HM, Belvedere O, Conway C, et al. Using next-generation sequencing for high resolution multiplex analysis of copy number variation from nanogram quantities of DNA from formalin-fixed paraffin-embedded specimens. *Nucleic Acids Res* 2010;38:e151.
- 38 Duncavage EJ, Magrini V, Becker N, et al. Hybrid capture and next-generation sequencing identify viral integration sites from formalin-fixed, paraffin-embedded tissue. *J Mol Diagn* 2011;13:325-33.
- 39 Chatterjee A, Leichter AL, Fan V, et al. A cross comparison of technologies for the detection of microRNAs in clinical FFPE samples of hepatoblastoma patients. *Sci Rep* 2015;5:10438.
- 40 Doleshal M, Magotra AA, Choudhury B, et al. Evaluation and validation of total RNA extraction methods for microRNA expression analyses in formalin-fixed, paraffin-embedded tissues. *J Mol Diagn* 2008;10(3):203-11.
- 41 Tumpey TM, Basler CF, Aguilar PV, et al. Characterization of the reconstructed 1918 Spanish influenza pandemic virus. *Science* 2005;310:77-80.
- 42 van den Brand JM, van Leeuwen M, Schapendonk CM, et al. Metagenomic analysis of the viral flora of pine marten and European badger feces. *J Virol* 2012;86:2360-65.
- 43 Bodewes R, van de Bildt MW, Schapendonk CM, et al. Identification and characterization of a novel adenovirus in the cloacal bursa of gulls. *Virology* 2013a;440:84-88.
- 44 Bodewes R, van Run P, Schürch AC, et al. Virus characterization and discovery in formalin-fixed paraffin-embedded tissues. *J Virol Meth* 2015;214:54-59.

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