

# Automation of Cell Staining

Technical Information Bulletin

## Automation of Cell Staining Using the Biomek 4000 Laboratory Automation Workstation

Amy N. Yoder, Sr. Development Scientist and Li Liu, Sr. Development Scientist, Beckman Coulter Life Sciences, Indianapolis, IN

### Abstract

Automation with the Biomek 4000 platform can improve the throughput and reproducibility of conventional staining of adherent cells. Applicability is demonstrated using mouse embryonic stem cells (mESCs) to detect the expression of developmental antigens in fixed cells following differentiation to cardiomyocytes which is relevant to stem cell biology research. The cardiomyocyte model resembles approaches used for the differentiation of other cell types. Automation protocols using the Biomek 4000 Laboratory Automation Workstation were developed using PerFix-nc fixation and permeabilization reagents with a no-wash protocol that provide a well resolved and reproducible detection of intracellular markers in murine embryonic stem cells and differentiated cardiomyocytes.

### Materials and Methods

#### *mESC Culture and Differentiation*

Mouse embryonic stem (mES) cells (Invitrogen, Carlsbad, CA) were maintained in growth media containing leukemia inhibitory factor (LIF) and 15% knockout serum replacement (KSR). For differentiation, cells were cultured in 15% FBS without LIF in a 384-well round bottom polypropylene plate (Nunc, Roskilde, Denmark) in 40  $\mu$ L differentiation medium (various treatments for 0 to 5 days). Embryoid bodies formed by Day 5 were transferred to gelatin-coated 96-well plates in 100  $\mu$ L fresh media. After Day 6, a portion of the adherent cells showed visible contraction. Control mES cells were treated in parallel without differentiation factors. Cells were harvested by Trypsin and dispensed into 96-well plates for processing.

#### *Biomek 4000 Laboratory Automation Workstation Configuration*

Figure 1 shows images of the instruments used in this information bulletin. Figure 2 shows the configuration of the Biomek 4000 Laboratory Automation Workstation used for cell staining experiments.







Figure 1. Workflow of automated microtiter plate-based cell staining process.

**Figure 2.** Automated cell staining method on the Biomek 4000 Laboratory Automation Workstation. **A)** Screen shot of the Define Pattern step in the automation method. Users can freely choose the sample wells on a 96-well plate format. **B)** Screen shot showing the Deck Layout. Tools, labware and tips are located as indicated. **C)** Screen shot of User Interface, where users can customize all the parameters, including reagent volumes, washes, incubation times and antibody information, etc. **D)** Screen shot showing the reagent volumes based on the parameters from the User Interface. **E)** Screen shot showing the custom defined antibody transfer step, where users can define the antibody locations manually.



File Edit Project Instrument Execution Options Help

New Open Close Save Preview Print Cut Copy Paste Undo Redo Clear Run Pause Stop Sim

Serial Dilution Aspirate Dispense Mix New Tip Unload Tip Loop Cleanup Move Pod Run Method Worklist Wait In Time Let If Script

Start  
Define Samples  
User Interface  
**Reservoir Volumes**  
Instrument Setup  
Wash Procedures  
Remove Media  
Fixation  
Permeabilization  
Blocking  
Staining  
Final Wash  
Final Processing  
Finish

### Cell Staining Reagent Volume for 36 Samples

8-Column Reagent Reservoir:

R1:	3120 µL	Fixation
R2:	3300 µL	Permeabilization
R3:	3300 µL	Blocking
R4:	3300 µL	Final Wash
R5:	0 µL	Mounting Media
R6:	Empty	Empty
R7:	Empty	Empty
R8:	Empty	Empty

8-Column Antibody Reservoir:

A1:	3300 µL	Ab1
A2:	3300 µL	Ab2
A3:	3300 µL	Ab3
A4:	3300 µL	Ab4
A5:	0 µL	Ab5
A6:	0 µL	Ab6
A7:	0 µL	Ab7
A8:	0 µL	Ab8

Full Modular Reservoir:

F1:	41000 µL	PBS
-----	----------	-----

D

File Edit Project Instrument Execution Options Help

New Open Close Save Preview Print Cut Copy Paste Undo Redo Clear Run Pause Stop Sim

Serial Dilution Aspirate Dispense Mix New Tip Unload Tip Loop Cleanup Move Pod Run Method Worklist Wait In Time Let If Script

Start  
Define Samples  
User Interface  
Reservoir Volumes  
Instrument Setup  
Wash Procedures  
Remove Media  
Fixation  
Permeabilization  
Blocking  
Staining  
If Staining  
Then  
If CustomSelect  
Then  
How to configure custom...  
Transfer = CustomVolume  
Transfer = CustomVolume  
Pause the whole system fo...  
Run removesupernatant  
If CustomWash  
End  
Else  
For i = 1 to =NumStain st...  
Scripted Let  
Transfer = Volur  
Pause the whole...

Used [Pod1] for transfer. Load Tool (AutoSelect)

Tip Handling  
Load 3300 µL tip and unload them when the transfer is done.  
Change tips between sources. Change tips between destinations.

Source: Antibody\_Tube  
Draw Well Contents from Antibody\_Tube using the Antibody Tube technique.

Destination: SamplePlate  
Dispense =CustomVolume of Tip Contents to SamplePlate using the Antibody technique.  
Click here to add a destination.

Stop when finished with Destinations  
Replicate each well 1 time  
Dispense up to 1 time per draw  
Aspirate at most 0 per transfer for repeated dispensing.

CellStainMethod CellStainProject\_Bome4000 STC: 6-49-24

E


### mES Cells and Cardiomyocytes Characterization

Differentiated cells were harvested on Day 8 using Accumax (Millipore, Billerica, MA). Then mES cells were mixed with either pooled differentiated cells or feeder cells (Invitrogen), and seeded into 96-well imaging plates and allowed 24 hours to attach. Using a Biomek 4000 Laboratory Automation Workstation with multi-probe tool, the mixed cells were fixed and permeabilized with PerFix-nc reagents (B10825, Beckman Coulter; Inc., Brea, CA). The mixture of differentiated cells and mES cells were stained with myosin heavy chain-Alexa Fluor 488 (clone MF20, eBioScience, San Diego, CA) and anti-Sox2-Alexa Fluor 647 (clone O30-678, BD Biosciences, San Jose, CA) (Figure 3). The mixture of mES cells and feeder cells were stained with either anti-Nanog-Alexa Fluor 488 (clone eBioMLC-51, eBioScience, San Diego, CA) or anti-Sox2-Alexa Fluor 647 (clone O30-678, BD Biosciences,

San Jose, CA). All conjugates were titered for optimal performance, and relevant isotype controls were used to control for non-specific staining. Cells were fixed by adding 5  $\mu$ L of PerFix-nc reagent 1 for 15 minutes, followed by permeabilization and staining with 50  $\mu$ L of PerFix-nc reagent 2 containing the antibody conjugates for 30 minutes. After the removal of the supernatant and replacement with 50  $\mu$ L reagent 3 for 5 minutes, cells were identified by nuclear staining with mounting medium which contains DAPI (Invitrogen).

### Imaging

All samples were analyzed on an ImageXpress system with MetaXpress software (Molecular Devices, Sunnyvale, CA).




**Figure 3.** Intracellular staining of stem cell marker Sox2 and myosin heavy chain using PerFix-nc reagents. PerFix-nc fixation and permeabilization reagents are optimized for no-wash detection of intracellular markers for flow cytometry, but also provide excellent sample preparation for imaging applications. The PerFix-nc kit was tested on undifferentiated and Day 8 differentiated mES cells for staining. Mixtures of mES cells and Day 8 differentiated mES cells (ratio 1:2) were detached and re-seeded into a 96-well imaging plate, then allowed 24 hours to plate out. PerFix-nc reagents and antibodies were added into the sample wells. All liquid handling steps were performed on a Biomek 4000 platform. Sixteen wells required less than an hour to prepare for imaging. Images were taken using the 10X objective on an ImageXpress high content imager and processed with MetaXpress software (Molecular Devices). Cells were stained with either isotype controls (B.a.) or with myosin heavy chain (green, A.b. and Sox2 (blue, A.c.). DAPI (red, A.d.) was used to mark the locations of nuclei. The group treated with cardiomyocyte-inducing reagents (B.c.) has more MF20 positive staining than the non-treated group (B.b.)

## Results

Expression of myosin heavy chain in differentiated cardiomyocytes (Figure 3) and stem cell markers (Figure 4) were clearly demonstrated in the cell preparation processed with the Biomek 4000 Laboratory Automation Workstation. Automated cell staining for 58 wells on the liquid handler required about

one hour without user intervention. The PerFix-nc sample preparation system provided good intracellular staining results with good reproducibility and accuracy without washing, greatly facilitating assay automation.



**Figure 4.** Mixtures of mES cells and feeder cells were seeded into a 96-well imaging plate, then allowed 24 hours to plate out before being treated with PerFix-nc reagents and stained with either anti-Nanog-Alexa Fluor 488 (A.b., green) or with anti-Sox2-Alexa Fluor 647 (A.c., blue). DAPI (A.d., red) was used to mark the locations of nuclei. The percentages of positive staining (B) were quantified using the MetaXpress software. The standard deviation coefficient of variation are the average values of 58 samples.

## Summary

This work demonstrates that cell staining sample preparation workflows can be automated using standard components on our new Biomek 4000 Laboratory Automation Workstation. Automation can achieve preparation timesavings with large

numbers of samples while maintaining equivalent results and precision compared with manual processing. As a result, it is possible to perform large cell staining studies in an expandable manner with walk-away capability.



\* The Biomek 4000 Workstation is available with an enclosure that is in development.  
All trademarks are the property of their respective owners.  
Beckman Coulter, the stylized logo and Biomek are trademarks of Beckman Coulter, Inc. and are registered with the USPTO.

For Beckman Coulter's worldwide office locations and phone numbers, please visit [www.beckmancoulter.com/contact](http://www.beckmancoulter.com/contact)